




μ FR[®] (AngioPlus Core)

--Angio-based FFR and PCI Planning

Pulse Medical



μFR analysis process

Angiography

Image selection: Recommended projection, clear angiography, less overlapping, less foreshortening

Flow

- Check flow velocity
- Select key frame
- Select ROI (Region Of Interest)
- Check vessel contour

2D μFR

- Branch detection
- Check reference lumen

2nd View

- Select key frame
- Select ROI
- Check vessel contour
- Corresponding marks

3D μFR

- Check 3D reconstruction morphology
- Check pressure pull back curve
- Check diameter change curve
- Check reference lumen

Export 2D analysis report

Export 3D analysis report

SOP of Coronary Angiographic Image Acquisition

1. Acquiring the angiographic image using recommended projection

Acquire **at least 1 angiographic image** runs according to the **recommended angiographic projections**. Adjust angiographic projections if the recommended projection has a suboptimal quality or vessel overlaps at the lesion segment.

2. Uniformly injection, fill the entire vessel

keep catheter to have better co-axiality with respect to vessel and maintain the uniformly contrast injection speed (~4 ml/s) until the distal vessel is filled with contrast and keep it to include the ED phase.

3. Nitroglycerin required

Inject nitroglycerin through the imaging/guiding catheter to **eliminate arterial spasm** prior to angiography.

4. Angiographic frame rate ≥ 15 frame/s, start recording before contrast injection, ≥ 1 cardiac cycle

Set angiographic imaging frame rate **≥ 15 frame/sec**. Then prefill the catheter with contrast medium. After complete wash out of the contrast medium, Start recording angiographic images at least **1 second prior to contrast injection**. The angiogram should include **≥ 1 cardiac cycle**, preferably 2 cycles.

5. Table moving

Don't move the table, especially the height.



2D Analysis - Six key principles for qualified μ FR analysis

Image Selection

- Recommended projection, clear angiography, less overlapping, less foreshortening

Key Frame

- Prioritize the frame with the maximal stenosis exposure
- Select the frame with the clearest vessel segment

Vessel Segmentation

- sufficient length ($\geq 20\text{mm}$), start point and end point in normal segment, moderate correction of contours

Flow Velocity

- Make sure the correct tracing of contrast/blood flow

Bifurcation







- Double check the detection and reference index position of bifurcations

Reference Lumen

- Ensure PN and DN in the positions with healthy lumen



2D Analysis - Image selection




Note: Preferably choose optimal angiography runs (especially in lesion segment) with less vascular overlaps and less vascular shortening.


2D Analysis - Key Frame

Degree of importance : $A > B = C > D$


A Maximal stenosis exposure




B Lower foreshortening



C Entire vessel clearly visualized



D Less overlapping






2D Analysis - Vessel Segmentation

Length: $\geq 20\text{mm}$ and includes all lesions

Start point: LM (LCA), beginning of the catheter (RCA) --(unless it is unclear and foreshortening or overlapping)

End point: normal segment with a marked point, exclude distortions



Vascular contours: good contrast- little correction; bad contrast- moderate correction



2D Analysis - Flow velocity

Check if the flow velocity curve, and adjust the range if necessary

- The grey box should cover all the frames, from the contrast entering the vessel, to the contrast flowing to the distal of vessel



Contains the period that contrast no longer flow further in the vessel (*platform period*, centerline length no longer increases)

correct flow velocity

2D Analysis - Flow velocity


If the centerline cannot be accurately detected but you still want to analyze, then the fixed flow model can be applied. But the μ FR results obtained from the fixed flow model should be regarded as for reference only.





2D Analysis - Bifurcation

Detection: identified all bifurcations (diameter $\geq$ 1.0mm): manually add missing branches, delete excess branches. Put the endpoint of branch after the last stenosis.


Reference Diameter: reference index should be in healthy segment.





Missed the big side-branch



Wrong endpoint





Wrong reference index



2D Analysis - Reference Lumen

PN -proximal normal **DN** -distal normal

Priority use of the automatic method



Fit well: reference lumen of the normal vessel segment fits the actual lumen (**Priority guarantee the normal segment adjacent to culprit lesion**)

PS: orthogonal reconstruction vessel will expand than normal segment


2D Analysis - Reference Lumen

PN -proximal normal **DN** -distal normal

Priority use of the automatic method

Mu (μ) values

1. Make sure all the large side branches are recognized correctly.
2. Check Mu (μ) values are within the normal range **0.7 - 1.3**.



3D Analysis - Four key principles for qualified μ FR analysis

Note: After the 2D analysis, please click "DONE" and then click "Save" to **save the 2D results**. Otherwise, once starting 3D analysis, it will be impossible to save the 2D results again.

Second Image Selection

- **Recommended angle difference $\geq 25^\circ$** , minimal difference $\geq 15^\circ$, clear angiography, less overlapping, less foreshortening

Key Frame

- Prioritize the frame with the maximal stenosis exposure
- Ensure the key frames of two selected image are in the same phase of the cardiac cycle

Coregistration

- **Start point and end point correspond to 1st image**
- Check the alignment of the two images
- Corresponding Marks

Check 3D Result

- Check 3D reconstruction morphology, reference lumen, pressure pull back curve

3D Analysis - Second Image Selection

Image Quality: optimal angiography (especially in lesion segment) with less overlapping, less foreshortening.


Second Image Angle: recommended angle difference $\geq 25^{\circ}$, minimal difference $\geq 15^{\circ}$.



3D Analysis - Key Frame

Ensure the key frames of two selected image are in the same phase of the cardiac cycle.



A. Maximal stenosis exposure > B. Lower foreshortening = C. Entire vessel clearly visualized > D. Less overlapping



3D Analysis - Coregistration


ROI: choose the same start point and end point with the 1st image.

Alignment: double check the alignment (especially at the most severe stenosis) and corresponding curve. If the contour is incorrectly identified, moderately adjust the contour of the second image.



3D Analysis - Coregistration

Corresponding Marks: If the alignment starts and ends correctly, but **the middle part** still cannot be matched, especially for the segment around the lesion, use the “corresponding marks” function to prioritize the alignment of the most severe stenosis.



3D Analysis - Check 3D Result

Check 3D reconstruction morphology, reference lumen, pressure pull back curve.

Double click to
change displayed
angiography image





博动 PULSE
MEDICAL

Demonstration



<http://en.pulse-imaging.com>



@PulseMedTech



Company/PulseMedical

